



FULL
PRODUCT
CATALOGUE




I. RADIOLOGY



Biliary Stents

Used as a palliative treatment for patients with biliary strictures.



01 Biliary Uncovered Stent

- Optimal radial force to help maintain patency.
- Unique cell construction provides good flexibility for placements in tortuous anatomies.

Art. No.	Diameter (mm)	Length (cm)	Delivery system
RV0804	8	4	7.5Fr / 50cm
RV0805		5	
RV0806		6	
RV0807		7	
RV0808		8	
RV0809		9	
RV0810		10	
RV0811		11	
RV0812		12	
RV1004	10	4	
RV1005		5	
RV1006		6	
RV1007		7	
RV1008		8	
RV1009		9	
RV1010		10	
RV1011		11	
RV1012		12	
RV1204	12	4	
RV1205		5	
RV1206		6	
RV1207		7	
RV1208		8	
RV1209		9	
RV1210		10	
RV1211		11	
RV1212		12	



02 Biliary Moving Cell Stent (MC Stent)

- Ideal for multiple stent-in-stent procedures.
- Unique moving cell design allows each cell of the stent to expand up to 10mm.

Art. No.	Diameter (mm)	Length (cm)	Delivery system
RVL0804	8	4	7.5Fr / 50cm
RVL0805		5	
RVL0806		6	
RVL0807		7	
RVL0808		8	
RVL0809		9	
RVL0810		10	
RVL0812		12	
RVL1004	10	4	
RVL1005		5	
RVL1006		6	
RVL1007		7	
RVL1008		8	
RVL1009		9	
RVL1010		10	
RVL1012		12	



03 Biliary Saessak Stent (Bud stent)

- Half covered Moving Cell stent.
- Uncovered Moving Cell portion enables multiple stent-in-stent procedures and allows better biliary flow from side branches.
- Covered PTFE portion lowers chances of tumor ingrowth.

Art. No.			Diameter (mm)	Length (cm)	Delivery system
40% covered	50% covered	60% covered			
RVLP0808-4	RVLP0808-5	RVLP0808-6	8	8	8.5Fr / 50cm
RVLP0810-4	RVLP0810-5	RVLP0810-6		10	
RVLP0812-4	RVLP0812-5	RVLP0812-6		12	
RVLP1008-4	RVLP1008-5	RVLP1008-6	10	8	
RVLP1010-4	RVLP1010-5	RVLP1010-6		10	
RVLP1012-4	RVLP1012-5	RVLP1012-6		12	



04 Biliary Covered Stent

- PTFE cover has been applied to prevent tumor ingrowth and create a smooth inner side to reduce sludge formation and food impaction.
- Wider heads were designed to reduce stent migration.

Art. No.		Diameter (mm)	Length (cm)	Stent head		Delivery system
Fully covered	Partially covered			Diameter (mm)	Length (cm)	
RVP0804HC	RVP0804HU	8	4	12	1 each	8.5Fr / 50cm
RVP0805HC	RVP0805HU		5			
RVP0806HC	RVP0806HU		6			
RVP0807HC	RVP0807HU		7			
RVP0808HC	RVP0808HU		8			
RVP0809HC	RVP0809HU		9			
RVP0810HC	RVP0810HU		10			
RVP0811HC	RVP0811HU		11			
RVP0812HC	RVP0812HU		12			
RVP1004HC	RVP1004HU	10	4	14		
RVP1005HC	RVP1005HU		5			
RVP1006HC	RVP1006HU		6			
RVP1007HC	RVP1007HU		7			
RVP1008HC	RVP1008HU		8			
RVP1009HC	RVP1009HU		9			
RVP1010HC	RVP1010HU		10			
RVP1011HC	RVP1011HU		11			
RVP1012HC	RVP1012HU		12			
RVP1204HC	RVP1204HU	12	4	16		
RVP1205HC	RVP1205HU		5			
RVP1206HC	RVP1206HU		6			
RVP1207HC	RVP1207HU		7			
RVP1208HC	RVP1208HU		8			
RVP1209HC	RVP1209HU		9			
RVP1210HC	RVP1210HU		10			
RVP1211HC	RVP1211HU		11			
RVP1212HC	RVP1212HU		12			



05 Biliary Spiral Cover Stent

- The PTFE cover in spiral configuration has been designed to allow better flexibility and biliary flow.
- Soft inner PTFE covered side reduces sludge formations.

Art. No.		Diameter (mm)	Length (cm)	Stent flange		Delivery system
Dumbbell type	T type			Diameter (mm)	Length (cm)	
RVN0804HC	RVN0804TC	8	4	12		8.5Fr / 50cm
RVN0805HC	RVN0805TC		5			
RVN0806HC	RVN0806TC		6			
RVN0807HC	RVN0807TC		7			
RVN0808HC	RVN0808TC		8			
RVN0809HC	RVN0809TC		9			
RVN0810HC	RVN0810TC		10			
RVN0811HC	RVN0811TC		11			
RVN0812HC	RVN0812TC		12			
RVN1004HC	RVN1004TC	10	4	14	Dumbbell type 1cm proximal & distal T type 1cm proximal	
RVN1005HC	RVN1005TC		5			
RVN1006HC	RVN1006TC		6			
RVN1007HC	RVN1007TC		7			
RVN1008HC	RVN1008TC		8			
RVN1009HC	RVN1009TC		9			
RVN1010HC	RVN1010TC		10			
RVN1011HC	RVN1011TC		11			
RVN1012HC	RVN1012TC		12			
RVN1204HC	RVN1204TC	12	4	16		
RVN1205HC	RVN1205TC		5			
RVN1206HC	RVN1206TC		6			
RVN1207HC	RVN1207TC		7			
RVN1208HC	RVN1208TC		8			
RVN1209HC	RVN1209TC		9			
RVN1210HC	RVN1210TC		10			
RVN1211HC	RVN1211TC		11			
RVN1212HC	RVN1212TC		12			



06 Biliary Jengu Stent (80cm long suture)

- Inner and outer PTFE cover help reduce tumor ingrowth and food impaction.
- 80cm long suture helps the stent to keep its place after stenting and be easily removed percutaneously.

Art. No.	Diameter (mm)	Length (cm)	Stent head		Delivery system
Fully covered			Diameter (mm)	Length (cm)	
RVP0804HC1-I	8	4	12	1 each	8.5Fr / 50cm
RVP0805HC1-I		5			
RVP0806HC1-I		6			
RVP0807HC1-I		7			
RVP0808HC1-I		8			
RVP0809HC1-I		9			
RVP0810HC1-I		10			
RVP0811HC1-I		11			
RVP0812HC1-I		12			
RVP1004HC1-I	10	4	14		
RVP1005HC1-I		5			
RVP1006HC1-I		6			
RVP1007HC1-I		7			
RVP1008HC1-I		8			
RVP1009HC1-I		9			
RVP1010HC1-I		10			
RVP1011HC1-I		11			
RVP1012HC1-I		12			
RVP1204HC1-I	12	4	16		
RVP1205HC1-I		5			
RVP1206HC1-I		6			
RVP1207HC1-I		7			
RVP1208HC1-I		8			
RVP1209HC1-I		9			
RVP1210HC1-I		10			
RVP1211HC1-I		11			
RVP1212HC1-I		12			

Contact us

sales@bcmglobal.co.kr

